






Sleep Review strives for accuracy in all data but cannot be held responsible for claims made by manufacturers. All apps may not be included. Email [sroy\[at\]medqor.com](mailto:sroy[at]medqor.com) to be considered for the next update.

Company		3B Medical Inc	Breas Medical Inc (formerly Human Design Medical)	Drive DeVilbiss Healthcare	Philips	ResMed
						
Mobile App		3B Luna QR	NiteLog	DeVilbiss SmartLink App	DreamMapper	myAir for Air10 by ResMed
Distribution Platforms	Google Play	X	X	X	X	X
	Apple App Store	X	X	X	X	X
	Amazon Appstore					
	BlackBerry World					
	Other		DataViewer website access	website access	website access	website access
PAP Device Compatibility		Luna, Luna II	Z1 CPAP, Z1 Auto	DV6 CPAPs (Bluetooth), DV5 (SmartCodes)	DreamStation, DreamStation Go, System One 50 & 60 Series	All Air10 devices
In-app Reporting	Apnea-Hypopnea Index	X	X	X	X	X
	Total Apnea-Hypopnea Count	X	X		X	
	Central Apnea Index	X		X		
	Total Central Apnea Count	X			X	
	Pressure Settings or Pressure Trends	X	X		X	
	Usage Hours	X	X	X	X	X
	Mask Leaks	X	X	X	X	X
Other	snore index, leak percentage, best 30			adherence score, % time in exhale puff (DV6x only), NOI (DV5x only), EPI (DV5x only), trends: usage AHI, leakage	hypopnea count, obstructed apnea count, comfort settings, mask disconnects, % of night in large leak	1-100 daily myAir score, broken down by usage hours, mask seal, mask on/off events, and apneas and hypopneas per hour
Available In-app Alerts/Notifications	Low Usage			X	X	X
	High Leak			X	X	X
	No Data Received			X	X	
	Resupply Soon			X		X
	Other				goal achievements, coaching messages, equipment cleaning reminders	30-day trends of usage hours, leak, AHI, mask on/off events, and myAir Score
In-app Patient Educational Resources		None, refers users to website	None	Educational videos	Videos, guides (images/text and PDFs). Content is defined by the equipment being used and the market where the account is registered.	Personalized videos and guides cover ResMed mask fit and cleaning (specific to patient's mask), equipment maintenance and comfort settings, and adjusting to therapy. Full FAQ.
HIPAA Compliant		X	X	X	X	X
Care Continuum Access		Healthcare team can access patient data at icodeconnect.com after registering at no cost.	Patient emails either data or screen shots to physician.	Evaluate data and reports through SmartLink Desktop and My SmartLink Cloud.	DreamMapper is connected to EncoreAnywhere/Care Orchestrator, which both provide EMR connectivity.	Access patients' device data through AirView, U-Sleep, or third-party data providers if integrated with ResMed.
Additional Information		None provided.	None provided.	The SmartLink App enables users to engage with IntelliPAP and IntelliPAP 2 CPAP devices. Users have an interactive role with their therapy and the SmartLink App enhances the user experience. The app captures quantitative data from the machine and qualitative data about how the patient is feeling on therapy.	DreamMapper is the top rated* and number-one downloaded** PAP therapy adherence app. *publicly available app store data **App Annie and Google Play Developer's Console data	Data from the Air10 is automatically uploaded daily. 87% of patients who used myAir and were monitored in AirView were PAP compliant by day 90 (CHEST 2017). Air10 device data populates in myAir and AirView automatically, once the user gives consent. Available on smartphones and online in 20+ countries.