



## EZ-ACCESS Proficiency Test

Business Name: \_\_\_\_\_ Name: \_\_\_\_\_  
Results: \_\_\_\_\_ Title: \_\_\_\_\_  
Date: \_\_\_\_\_

*Carefully read each question and circle the correct answer 80% required for completion*

EZ-ACCESS SOLUTION BASED PRODUCTS MAKE SAFE INDEPENDENT LIVING A POSSIBILITY

1. What Ez-Access Product allows to make patient transfers and activities easier for both patient and their caregiver. This Patient Transfer device pivots 360 degrees for complete coverage.

- a. SUITCASE Portable Suitcase Ramp
- b. ORBIT Patient Transfer Lift
- c. TILT Toilet Incline Lift



2. The ORBIT Patient Transfer Lift can be used as an alternative to Ceiling or Portable Track System for Patient Transfers?

- a. Yes
- b. No

3. What Ez-Access Product allows both safety and comfort as a lift chair without upholstery that offers safe and independent bathroom use by assisting in getting on and off the commode?

- a. Vertical Platform Lift "VPL"
- b. ORBIT Patient Transfer Lift
- c. TILT Toilet Incline Lift



4. What types of commodes does the TILT Toilet Incline Lift accommodate?

- a. Standard Toilet
- b. Elongated Toilet
- c. Offers both a Standard and Elongated Toilet System for both types of commodes

5. What Ez-Access Product allows for a smooth transition from the floor to the door to eliminate potential trip hazards. Designed for doorways, sliding glass doors, and raised landings.

- a. Pathway 3G Modular System
- b. TRANSITIONS Threshold Ramps
- c. GATEWAY Solid Surface Portable Ramps



6. How do you measure for an EZ-ACCESS TRANSITIONS Threshold Ramp?

- a. Distance in height in inches from the floor to the door
- b. Measurement from the front of the door to the back of the door
- c. The width of the door being measured

7. What is the most aesthetically pleasing TRANSITIONS Threshold Ramps and what colors are available? HINT: This is the most popular.

- a. Modular Entry Ramp in aluminum silver
- b. Angled Entry Mat in Black, Grey or Brown
- c. Angled Entry Plate in Black



8. What Ez-Access Product is the ideal for overcoming curbs, steps, low landings and loading mobility products into a vehicle? They are designed to setup quickly and folds for portability.

- a. GATEWAY Solid Surface Portable Ramps
- b. TRANSITIONS Threshold Ramps
- c. SUITCASE Portable Folding Ramps



9. What is a SUITCASE Top Lip Extension?

- a. It is a three-foot ramp extension for a portable ramp
- b. It extends the top transition plate to 9" allowing the ramp to clear the rear bumper when loading a mobility device into the back of an SUV or van
- c. It is a bumper that is added to the ramp to provide a soft landing in case of a fall



10. What Type of SUITCASE Portable Ramp allows for ultimate portability when longer ramp lengths are needed for easier vehicle storage?

- a. GATEWAY 3G Solid Surface Portable Ramps
- b. PATHWAY 3G Modular Access System Ramps
- c. SUITCASE TRIFOLD Advantage Series Ramps



11. What is the difference between Signature Series and Advantage Series Portable Ramps?

- a. Signature Series Ramps offer extruded slip resistance aluminum surface, and Advantage Series Ramps have applied slip resistance surfaces
- b. Signature Series Ramps are no different from Advantage Series Ramps
- c. Signature Series Ramps have applied slip resistance surfaces, and Advantage Series Ramps offer extruded slip resistance aluminum surface

12. What is Warranty for all styles of Aluminum Ramps and Threshold Ramps

- a. 3-year warranty
- b. Limited Lifetime Warranty
- c. 25-year warranty

13. Where are Ez-Access Aluminum Ramps and Threshold Ramp Manufactured

- a. China
- b. United States – American Made
- c. Mexico

14. What Ez-Access product allows safe and seamless access for both standard and unique entryways? This Ramping is designed for residential use when accessibility needs extend beyond the scope of a Portable ramping system. This system is semi-permanent and can be customized as needed.

- a. GATEWAY 3G Solid Surface Portable Ramp
- b. PATHWAY 3G Modular Access Ramping System
- c. PASSPORT Vertical Platform Lift



15. What is the most important question when discussing PATHWAY 3G Modular Ramping Systems?
  - a. What is your budget for your Home Accessibility Needs
  - b. Does the Customer want to use the Ramping System on their own, or will they have assistance of a caregiver
  - c. What is the total Rise from the ground to the door to measure the length of ramping needed
  
16. What does the American Disability Act (ADA) recommend in length of ramp to use the ramp for self-use?
  - a. 1:1 ratio (for every inch of rise, one foot of ramp is recommended)
  - b. 1:2 ratio (for every inch of rise, two feet of ramp is recommended)
  - c. 1:10 ratio (for every inch of rise, ten feet of ramp is recommended)
  
17. What type of customer requires us to use "ADA" guidelines when quoting ramp lengths?
  - a. Consumer
  - b. Businesses and State Funded Agencies
  - c. Social Workers
  
18. Does a consumer have to purchase a ramp that meets "ADA" guidelines?
  - a. Yes, they need to purchase for the safest outcome
  - b. No, a consumer can order whatever they want, the dealer just advises what is the safest outcome for the patient such as "ADA" guidelines
  
19. What is the highest Rise recommended for a patient that has the assistance of a caregiver?
  - a. 1:1 ration (for every inch of rise, one foot of ramp is recommended)
  - b. 1:2 ratio (for every inch of rise, two feet of ramp is recommended)
  - c. 1:10 ratio (for every inch of rise, ten feet of ramp is recommended)

20. What Ez-Access product is the ideal solution when surrounding space is limited, or higher than 30" of total rise are present for in home for access needs?

- a. PATHWAY 3G Modular Ramping System
- b. PASSPORT Vertical Platform Lift
- c. GATEWAY 3G Solid Surface Portable Ramp



21. Who is the ideal customer for a Modular Ramping System or Vertical Platform Lift application?

- a. Patient with a rollator
- b. Patient with a Power Wheelchair
- c. Patient with a Cane

22. How is a PASSPORT Vertical Platform Lift Powered?

- a. AC Powered (Direct Plug In)
- b. DC Powered (24V Batteries for backup purposes with Corded Trickle Charger)
- c. Both

23. How many cycles can a PASSPORT Vertical Platform Lift run on just the 24V Batteries?

- a. 10 complete cycles in case of power failure
- b. 1 complete cycle in case of power failure
- c. 30 complete cycles in case of a power failure

24. What is weight Capacity of PASSPORT Vertical Platform Lift?

- a. 500 Lbs
- b. 750 Lbs
- c. 1000 Lbs

25. How Long has Ez-Access been in business providing access to life beyond barriers?

- a. 10+ Years
- b. 20+ Years
- c. 30+ Years